

Mounting Stand

MS-7 Telescoping Stand



Mounting stands for gooseneck fiber bundles are used to help you position fiber bundles to get light in the spot you need. Mounting stands give you the option to keep your illuminator further away from your work area.

The MS-7 Telescoping Stand is typically used with the EEG3736 fiber optic light guide. Use of two bundles for opposed angles serves to diminish shadows, reduce optical axis glare, and provides a constant field of illumination.

Individuals can control intensities by varying the distance of the fiber optic light output to the object plane.